

Anti-Card14 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bxl6

Item Number: CM0022220

Introduction

CARD14 rabbit polyclonal antibody validated for Western blot, immunohistochemistry (paraffin), and ELISA. Detects human, mouse, and rat CARD14. Ideal for studying NF-κB inflammatory signaling and related pathways.

[Learn More](#)

Parameter	Value
Product Name	CARD14 Rabbit pAb
Remarks/Alias	PRP; PSS1; BIMP2; CARMA2; PSORS2; CARD14
Species	Human
GeneID (Human)	79092
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-500 of human CARD14 (NP_001244899.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXL6
Protein Weight	113kDa
Shipping	Ice bag