


Anti-Card19 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Lw7

Item Number: CM0015367



Introduction

High-quality anti-CARD19 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and mouse CARD19, a BCL10-interacting protein involved in NF-kappa-B inhibition. Suitable for cell signaling and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	CARD19 Rabbit pAb
Remarks/Alias	BinCARD; C9orf89; CARD19
Species	Human
Gene ID (Human)	84270
Gene ID	84270
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CARD19 (NP_115686.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96LW7
Protein Weight	26kDa
Shipping	Ice bag