

Anti-Card8 Polyclonal Antibody For Wb, Elisa - Q9Y2G2

Item Number: CM0014217



Introduction

CARD8 Rabbit pAb for WB and ELISA applications. Detects human caspase recruitment domain-containing protein 8, a key inflammasome sensor. Suitable for pyroptosis, innate immunity, and HIV-1 latency studies.

[Learn More](#)

Parameter	Value
Product Name	CARD8 Rabbit pAb
Remarks/Alias	NDPP; DACAR; DAKAR; NDPP1; TUCAN; CARDINAL; CARD8
Species	Human
GeneID (Human)	22900
GeneID	22900
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CARD8 (NP_001171829.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y2G2
Protein Weight	61 kDa
Shipping	Ice bag