


Anti-Cardiac Troponin I (Tnni3) Polyclonal Antibody For Wb, Ip, If-P, Elisa - P19429

Item Number: CM0013916



Introduction

High-quality rabbit polyclonal antibody targeting human Cardiac Troponin I (TNNI3). Suitable for WB, IP, IF-P, and ELISA with cross-reactivity to mouse and rat. Ideal for cardiac muscle research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Cardiac Troponin I (TNNI3) Rabbit pAb
Remarks/Alias	CMH7; RCM1; cTnI; CMD2A; TNNC1; CMD1FF; Cardiac Troponin I (TNNI3)
Species	Human
GeneID (Human)	7137
GeneID	7137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human Cardiac Troponin I (TNNI3) (NP_000354.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P19429
Protein Weight	24kDa
Shipping	Ice bag