


Anti-Cardiac Troponin I (Tnni3) Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P19429

Item Number: CM0003163



Introduction

High-quality rabbit monoclonal antibody against human Cardiac Troponin I (TNNI3), validated for WB, IF-P, ELISA; cross-reacts with mouse and rat. Ideal for cardiac research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Cardiac Troponin I (TNNI3) Rabbit mAb
Remarks/Alias	CMH7; RCM1; cTnI; CMD2A; TNNC1; CMD1FF; Cardiac Troponin I (TNNI3)
Species	Human
GeneID (Human)	7137
GeneID	7137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human Cardiac Troponin I (TNNI3) (NP_000354.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P19429
Protein Weight	24kDa
Shipping	Ice bag