

Anti-Cardiac Troponin T (Tnnt2) Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P45379

Item Number: CM0005231



Introduction

Rabbit monoclonal antibody targeting human cardiac troponin T (TNNT2), validated for IF-P, IHC-P, ELISA. Detects human, mouse, and rat. Ideal for cardiac research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Cardiac troponin T (TNNT2) Rabbit mAb
Remarks/Alias	CMH2; RCM3; TnTC; cTnT; CMD1D; CMPD2; LVNC6; Cardiac troponin T (TNNT2)
Species	Human
Gene ID (Human)	7139
Gene ID	7139
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45379
Protein Weight	31-36 kDa
Shipping	Ice bag