

Anti-Cardiac Troponin T/Tnnt2 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P45379

Item Number: CM0011557



Introduction

Rabbit monoclonal antibody against human cardiac troponin T (TNNT2), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cardiac research.

[Learn More](#)

Parameter	Value
Product Name	Cardiac Troponin T/TNNT2 Rabbit mAb
Remarks/Alias	CMH2, RCM3, TnTC, cTnT, CMD1D, CMPD2, LVNC6
Species	Human
Gene ID (Human)	7139
Gene ID	7139
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-288 of human cTNT (NP_001001430.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45379
Protein Weight	31-36 kDa
Shipping	Ice bag