

Anti-Cars Polyclonal Antibody (Rabbit) For Wb, If/Icc, Elisa - P49589

Item Number: CM0029807

Introduction

CARS Rabbit pAb is a polyclonal antibody for WB, IF/ICC, ELISA applications. It detects human and mouse CARS, an 85kDa cytoplasmic cysteine--tRNA ligase essential for protein biosynthesis. Ideal for translation research.

[Learn More](#)

Parameter	Value
Product Name	CARS Rabbit pAb
Remarks/Alias	CARS; MDBH; CYSRS; MCDDBH; MGC:11246
Species	Human
Gene ID (Human)	833
Gene ID	833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 469-748 of human CARS (NP_001742.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49589
Protein Weight	85kDa
Shipping	Ice bag