

Anti-Cart Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16568

Item Number: CM0028218



Introduction

Rabbit polyclonal antibody targeting human cocaine- and amphetamine-regulated transcript (CART) protein. Validated for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for obesity, feeding, and neuropeptide research.

[Learn More](#)

Parameter	Value
Product Name	CART Rabbit pAb
Remarks/Alias	CART
Species	Human
Gene ID	9607
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human CART (NP_004282.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16568
Protein Weight	13kDa
Shipping	Ice bag