

Anti-Cartpt Polyclonal Antibody For Ihc-P, Elisa - Q16568

Item Number: CM0021328



Introduction

High-quality rabbit polyclonal antibody against human CARTPT (CART) for IHC-P and ELISA applications. Cross-reacts with mouse and rat. Target protein weight 13kDa. SWISS: Q16568.

[Learn More](#)

Parameter	Value
Product Name	CARTPT Rabbit pAb
Remarks/Alias	CART; CARTPT
Species	Human
GeneID (Human)	9607
GeneID	9607
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CARTPT (NP_004282.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16568
Protein Weight	13kDa
Shipping	Ice bag