


Anti-Cas9 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q99Zw2

Item Number: CM0010898



Introduction

Rabbit monoclonal Cas9 antibody validated for WB, IF/ICC, ELISA. Recognizes Streptococcus pyogenes Cas9 (Q99ZW2), predicted to cross-react with human and other species. Ideal for CRISPR research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Cas9 Rabbit mAb
Remarks/Alias	CAS9
Species	Streptococcus pyogenes
GeneID	901176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-730 of Streptococcus pyogenes M1 GAS CAS9 (NP_269215.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Other (Wide Range Predicted)
SWISS	Q99ZW2
Protein Weight	158kDa
Shipping	Ice bag