

Anti-Casein Kinase 1 Alpha (Csnk1A1) Monoclonal Antibody For Wb, Ihc-P, Elisa - P48729

Item Number: CM0011463



Introduction

High-quality rabbit monoclonal antibody targeting human Casein Kinase 1 alpha (CSNK1A1). Validated for WB, IHC-P, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for Wnt signaling, cell cycle, and kinase research applications.

[Learn More](#)

Parameter	Specification
Product Name	Casein Kinase 1 alpha (CSNK1A1) Rabbit mAb
Remarks/Alias	CK1; CK1a; CKIa; HLCDGP1; PRO2975; HEL-S-77p; Casein Kinase 1 alpha (CSNK1A1)
Species	Human
GeneID (Human)	1452
GeneID	1452
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 238-337 of human Casein Kinase 1 alpha (CSNK1A1) (P48729).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48729
Protein Weight	39kDa
Shipping	Ice bag