

## Anti-Casein Kinase 1 Alpha (Csnk1A1) Polyclonal Antibody For Wb - P48729

Item Number: CM0019836



### Introduction

Rabbit polyclonal antibody against human CSNK1A1, suitable for WB, IHC-P, IF/ICC, ELISA. Validated in human, mouse, and rat samples for studying Wnt signaling and kinase activity.

[Learn More](#)

| Parameter        | Value                                                                                                                           |
|------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Casein Kinase 1 alpha (CSNK1A1) Rabbit pAb                                                                                      |
| Remarks/Alias    | CK1, CK1a, CKIa, HLCDGP1, PRO2975, HEL-S-77p                                                                                    |
| Species          | Human                                                                                                                           |
| Gene ID (Human)  | 1452                                                                                                                            |
| Gene ID          | 1452                                                                                                                            |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 238-337 of human Casein Kinase 1 alpha (CSNK1A1) (NP_001883.4) |
| Source           | Rabbit                                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                                           |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                               |
| SWISS            | P48729                                                                                                                          |
| Protein Weight   | 39kDa                                                                                                                           |
| Shipping         | Ice bag                                                                                                                         |