

Anti-Casein Kinase 2 Alpha (Csnk2A1) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P68400

Item Number: CM0023583



Introduction

High-quality rabbit polyclonal antibody against human Casein Kinase 2 alpha (CSNK2A1), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studying cell cycle, apoptosis, transcription, and DNA damage response.

[Learn More](#)

Parameter	Specification
Product Name	Casein Kinase 2 alpha (CSNK2A1) Rabbit pAb
Remarks/Alias	CKII; Cka1; Cka2; CK2A1; OCNDS; Casein Kinase 2 alpha (CSNK2A1)
Species	Human
Gene ID (Human)	1457
Gene ID	1457
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-391 of human Casein Kinase 2 alpha (CSNK2A1) (NP_001886.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68400
Protein Weight	45kDa
Shipping	Ice bag