

Anti-Casein Kinase 2 Beta (Csnk2B) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P67870

Item Number: CM0008266



Introduction

Validated recombinant rabbit monoclonal antibody against human Casein Kinase 2 beta (CSNK2B), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying kinase signaling, EBV infection, and neurodevelopmental disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Casein Kinase 2 beta (CSNK2B) Rabbit mAb
Remarks/Alias	G5A; CK2B; CK2N; Ckb1; Ckb2; CSK2B; POBINDS; Casein Kinase 2 beta (CSNK2B)
Species	Human
GeneID (Human)	1460
GeneID	1460
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 101-215 of human Casein Kinase 2 beta (CSNK2B) (CSNK2B) (P67870).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67870
Protein Weight	25kDa
Shipping	Ice bag