

Anti-Casein Kinase 2 Beta (Csnk2B) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P67870

Item Number: CM0022947



Introduction

Rabbit polyclonal antibody against human Casein Kinase 2 beta (CSNK2B) validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying CK2 signaling, Wnt pathway, and EBV infection research.

[Learn More](#)

Parameter	Value
Product Name	Casein Kinase 2 beta (CSNK2B) Rabbit pAb
Remarks/Alias	G5A; CK2B; CK2N; Ckb1; Ckb2; CSK2B; POBINDS
Species	Human
Gene ID (Human)	1460
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human Casein Kinase 2 beta (CSNK2B) (NP_001311.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67870
Protein Weight	25kDa
Shipping	Ice bag