

Anti-Cask Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O14936

Item Number: CM0021181



Introduction

CASK rabbit polyclonal antibody (pAb) raised against recombinant human CASK (aa 1-300). For WB, IF/ICC, ELISA. Ideal for synaptic protein scaffolding and kinase signaling research. Human reactive. UniProt: O14936.

[Learn More](#)

Parameter	Value
Product Name	CASK Rabbit pAb
Remarks/Alias	CMG; FGS4; LIN2; TNRC8; hCASK; CAGH39; CAMGUK; MICPCH; MRXSNA; CASK
Species	Human
GeneID (Human)	8573
GeneID	8573
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O14936
Protein Weight	105kDa
Shipping	Ice bag