


## Anti-Caspase-1 Monoclonal Antibody For Wb, Elisa - P29466

Item Number: CM0008295



### Introduction

Rabbit monoclonal antibody against human Caspase-1 (ICE, p45), validated for Western blot and ELISA. Targets the inflammatory caspase involved in IL-1 $\beta$ /IL-18 processing and pyroptosis. Suitable for immunology and cell death research.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Caspase-1 Rabbit mAb                                                                                                   |
| Remarks/Alias    | ICE; P45; IL1BC                                                                                                        |
| Species          | Human                                                                                                                  |
| Gene ID (Human)  | 834                                                                                                                    |
| Gene ID          | 834                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 120-297 of human Caspase-1 (NP_150634.1) |
| Source           | Rabbit                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                  |
| Application      | WB, ELISA                                                                                                              |
| Cross Reactivity | Human                                                                                                                  |
| SWISS            | P29466                                                                                                                 |
| Protein Weight   | 45kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |