


Anti-Caspase-1 Polyclonal Antibody For Wb, If-P, Elisa - P29466

Item Number: CM0022380



Introduction

Caspase-1 (ICE) rabbit polyclonal antibody for WB, IF-P, ELISA. Detects human, mouse, rat. Ideal for inflammasome-mediated IL-1 β processing, pyroptosis, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Caspase-1 Rabbit pAb
Remarks/Alias	ICE; P45; IL1BC; CASP1
Species	Human
Gene ID (Human)	834
Gene ID	834
Immunogen	Synthetic peptide corresponding to amino acids 198-297 of human CASP1 (NP_150634.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29466
Protein Weight	45kDa
Shipping	Ice bag