

Anti-Caspase-1 Rabbit Monoclonal Antibody For Wb, Elisa - P29452

Item Number: CM0005974

Introduction

Caspase-1 Rabbit mAb (P29452) validated for WB, ELISA in mouse. Recombinant antibody targeting amino acids 119-296. Essential for studying inflammatory responses and pyroptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	Caspase-1 Rabbit mAb
Remarks/Alias	ICE; IL1bc
Species	Mouse
GeneID (Human)	834
GeneID	12362
Immunogen	Recombinant protein corresponding to a sequence within amino acids 119-296 of mouse CASP1 (NP_033937.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P29452
Protein Weight	46 kDa
Shipping	Ice bag