

Anti-Caspase-1 Rabbit Polyclonal Antibody For Wb And Elisa - P29452

Item Number: CM0021772

Introduction

Rabbit polyclonal antibody targeting mouse Caspase-1 (ICE). Validated for Western blot and ELISA. Detects a ~46 kDa protein highly expressed in spleen and lung. Suitable for inflammation and pyroptosis research.

[Learn More](#)

Parameter	Specification
Product Name	Caspase-1 Rabbit pAb
Remarks/Alias	ICE; IL1bc; CASP1
Species	Mouse
GeneID (human)	834
GeneID	12362
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of mouse CASP1 (NP_033937.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P29452
Protein Weight	46kDa
Shipping	Ice bag