


Anti-Caspase-1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P29466

Item Number: CM0022083



Introduction

Rabbit polyclonal antibody against human Caspase-1 (ICE, IL1BC). Validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Suitable for research on inflammation, pyroptosis, and IL-1 β processing.

[Learn More](#)

Parameter	Value
Product Name	Caspase-1 Rabbit pAb
Remarks/Alias	ICE; P45; IL1BC
Species	Human
GeneID (Human)	834
GeneID	834
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-297 of human Caspase-1 (NP_150634.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29466
Protein Weight	45kDa
Shipping	Ice bag