


Anti-Caspase-1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P29466

Item Number: CM0016000



Introduction

Rabbit polyclonal antibody targeting human Caspase-1, validated for WB, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Ideal for research on inflammation and pyroptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	Caspase-1 Rabbit pAb
Remarks/Alias	ICE; P45; IL1BC; Caspase-1
Species	Human
Gene ID	834
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 201-300 of human Caspase 1 (NP_150634.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29466
Protein Weight	45kDa
Shipping	Ice bag