

Anti-Caspase-12 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Uxs9

Item Number: CM0018667

Introduction

Caspase-12 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Negative regulator of inflammatory responses; ideal for innate immunity and NF-kB pathway studies. - Q6UXS9.

[Learn More](#)

Parameter	Value
Product Name	Caspase-12 Rabbit pAb
Remarks/Alias	CASP-12; CASP12P1; Caspase-12
Species	Human
GeneID (human)	100506742
GeneID	100506742
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Caspase-12 (NP_001177945.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UXS9
Protein Weight	39kDa
Shipping	Ice bag