

Anti-Caspase-14 Monoclonal Antibody For Wb, Elisa - P31944

Item Number: CM0004814

Introduction

Validated for Western blot and ELISA, this Caspase-14 Rabbit monoclonal antibody detects human, mouse, and rat targets, supporting research into epidermal differentiation, cornification, and keratinocyte function.

[Learn More](#)

Parameter	Value
Product Name	Caspase-14 Rabbit mAb
Remarks/Alias	ARCI12; Caspase-14
Species	Human
GeneID (human)	23581
GeneID	23581
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human Caspase-14 (P31944).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31944
Protein Weight	28kDa
Shipping	Ice bag