

Anti-Caspase-14 Polyclonal Antibody For Wb, Elisa - P31944

Item Number: CM0027597

Introduction

Caspase-14 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for epidermal differentiation and keratinocyte apoptosis studies. 28kDa protein weight.

[Learn More](#)

Parameter	Value
Product Name	Caspase-14 Rabbit pAb
Remarks/Alias	ARCI12; Caspase-14
Species	Human
Gene ID (Human)	23581
Gene ID	23581
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human Caspase-14 (NP_036246.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31944
Protein Weight	28kDa
Shipping	Ice bag