


Anti-Caspase-2 Monoclonal Antibody For Wb, Elisa - P42575

Item Number: CM0004574



Introduction

Caspase-2 Rabbit Monoclonal Antibody validated for WB and ELISA in human samples. Detects 51kDa proenzyme. Key regulator of apoptosis via PIDDosome complex. Research use only.

[Learn More](#)

Parameter	Value
Product Name	Caspase-2 Rabbit mAb
Remarks/Alias	ICH1; NEDD2; CASP-2; NEDD-2; PPP1R57; Caspase-2
Species	Human
GeneID (Human)	835
GeneID	835
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 334-452 of human Caspase-2 (P29594).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42575
Protein Weight	51kDa
Shipping	Ice bag