


Anti-Caspase-2 Rabbit Polyclonal Antibody For Wb And Elisa - P42575

Item Number: CM0028584



Introduction

Caspase-2 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA applications. It reacts with human, mouse, and rat samples and detects the 51 kDa caspase-2 protein. Suitable for apoptosis and PIDDosome pathway research.

[Learn More](#)

Parameter	Value
Product Name	Caspase-2 Rabbit pAb
Remarks/Alias	ICH1; NEDD2; CASP-2; NEDD-2; PPP1R57; Caspase-2
Species	Human
Gene ID (Human)	835
Gene ID	835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 153-452 of human Caspase-2 (NP_116764.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42575
Protein Weight	51kDa
Shipping	Ice bag