


Anti-Caspase-3 Monoclonal Antibody For Wb, Elisa - P42574

Item Number: CM0003942



Introduction

KO-validated rabbit monoclonal antibody targeting Caspase-3, a key executioner caspase in apoptosis. Suitable for Western blot and ELISA in human and mouse samples. Ideal for apoptosis and cell death research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Caspase-3 Rabbit mAb
Remarks/Alias	CPP32; SCA-1; CPP32B; -3
Species	Human
Gene ID (Human)	836
Gene ID	836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-175 of human Caspase-3 (P42574)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42574
Protein Weight	32kDa
Shipping	Ice bag