


Anti-Caspase-3 Monoclonal Antibody For Wb, Elisa - P42574

Item Number: CM0003752



Introduction

Caspase-3 mouse monoclonal antibody validated for WB and ELISA, specific to human and mouse caspase-3. Ideal for apoptosis pathway analysis, detection of cleaved caspase-3, and programmed cell death research.

[Learn More](#)

Parameter	Value
Product Name	Caspase-3 Mouse mAb
Remarks/Alias	CPP32; SCA-1; CPP32B; Caspase-3
Gene ID (Human)	836
Gene ID	836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-160 of human Caspase-3 (P42574).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42574
Protein Weight	32kDa
Shipping	Ice bag