


Anti-Caspase-3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42574

Item Number: CM0024244



Introduction

Rabbit polyclonal antibody against human Caspase-3, suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for apoptosis and effector caspase studies.

[Learn More](#)

Parameter	Value
Product Name	Caspase-3 Rabbit pAb
Remarks/Alias	CPP32; SCA-1; CPP32B; -3
Species	Human
GeneID (Human)	836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-160 of human Caspase-3 (NP_004337.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42574
Protein Weight	32kDa
Shipping	Ice bag