


Anti-Caspase-3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42574

Item Number: CM0017609



Introduction

CamelBio offers rabbit polyclonal anti-Caspase-3 antibody (CPP32, SCA-1) for apoptosis research. Validated in WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Essential tool for studying programmed cell death.

[Learn More](#)

Parameter	Value
Product Name	Caspase-3 Rabbit pAb
Remarks/Alias	CPP32; SCA-1; CPP32B; Caspase-3
Species	Human
GeneID	836
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 6-105 of human Caspase-3 (NP_004337.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42574
Protein Weight	32kDa
Shipping	Ice bag