

Anti-Caspase-4 Rabbit Pab For Wb, Ihc-P, Elisa - P49662

Item Number: CM0024530



Introduction

Caspase-4 Rabbit pAb (TX, Mih1, ICH-2) validated for WB, IHC-P, ELISA in human. Targets inflammatory caspase involved in non-canonical inflammasome-mediated pyroptosis and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Caspase-4 Rabbit pAb
Remarks/Alias	TX, Mih1, ICH-2, Mih1/TX, ICEREL-II, ICE(rel)II, Caspase-4
Species	Human
Gene ID	837
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 81-270 of human Caspase-4 (NP_001216.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P49662
Protein Weight	43kDa
Shipping	Ice bag