


Anti-Caspase-6 Monoclonal Antibody For Wb, Elisa - P55212

Item Number: CM0005434



Introduction

High-quality anti-Caspase-6 rabbit monoclonal antibody validated for WB and ELISA. Detects human, mouse, rat caspase-6. Ideal for apoptosis, neurobiology, and innate immunity research. - P55212.

[Learn More](#)

Parameter	Value
Product Name	Caspase-6 Rabbit mAb
Remarks/Alias	MCH2; CSP-6; caspase-6; Caspase-6
Species	Human
GeneID (Human)	839
GeneID	839
Immunogen	Synthetic Peptide: A synthesized peptide derived from human Caspase-6.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55212
Protein Weight	33kDa
Shipping	Ice bag