


Anti-Caspase-6 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P55212

Item Number: CM0007475



Introduction

Rabbit monoclonal Caspase-6 antibody for WB, IHC-P, ELISA detects 33 kDa protein across human, mouse, and rat. Suitable for apoptosis and virology research applications, including PANoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	Caspase-6 Rabbit mAb
Remarks/Alias	MCH2; CSP-6; caspase-6; Caspase-6
Species	Human
GeneID (Human)	839
GeneID	839
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human Caspase-6 (P55212).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55212
Protein Weight	33kDa
Shipping	Ice bag