


Anti-Caspase-6 Rabbit Monoclonal Antibody For Wb - P55212

Item Number: CM0002744



Introduction

Caspase-6 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human and rat. Ideal for apoptosis, cell death, and neurodegeneration research. - P55212

[Learn More](#)

Parameter	Value
Product Name	Caspase-6 Rabbit mAb
Remarks/Alias	MCH2; CSP-6; caspase-6; Caspase-6
Species	Human
GeneID (human)	839
GeneID	839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-162 of human Caspase-6 (P55212).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P55212
Protein Weight	33 kDa
Shipping	Ice bag