


Anti-Caspase-6 Rabbit Polyclonal Antibody For Wb, Elisa - P55212

Item Number: CM0020029



Introduction

Rabbit polyclonal antibody against human Caspase-6 (SWISS P55212). Validated for WB and ELISA with cross-reactivity to human and mouse. Essential for apoptosis, axonal degeneration, and innate immunity research.

[Learn More](#)

Parameter	Specification
Product Name	Caspase-6 Rabbit pAb
Remarks/Alias	MCH2, CSP-6, caspase-6, Caspase-6
Species	Human
GeneID (Human)	839
Gene ID	839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 194-293 of human Caspase-6 (NP_001217.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P55212
Protein Weight	33kDa
Shipping	Ice bag