


Anti-Caspase-7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P55210

Item Number: CM0016625



Introduction

Caspase-7 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Specific to the human Caspase-7 protein for apoptosis and cell death research.

[Learn More](#)

Parameter	Value
Product Name	Caspase-7 Rabbit pAb
Remarks/Alias	MCH3; CMH-1; LICE2; CASP-7; ICE-LAP3; Caspase-7
Species	Human
Gene ID (Human)	840
Gene ID	840
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human Caspase-7 (NP_001218.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55210
Protein Weight	34kDa
Shipping	Ice bag