

Anti-Caspase-9 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P55211

Item Number: CM0019252



Introduction

High-quality Caspase-9 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat; detects 46 kDa protein. Ideal for apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	Caspase-9 Rabbit pAb
Remarks/Alias	MCH6; APAF3; APAF-3; PPP1R56; ICE-LAP6; Caspase-9
Species	Human
Gene ID (Human)	842
Gene ID	842
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 139-416 of Caspase-9 (NP_001220.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55211
Protein Weight	46kDa
Shipping	Ice bag