

Anti-Caspr1 Rabbit Monoclonal Antibody For Wb And Elisa - P78357

Item Number: CM0011302

Introduction

High-quality anti-Caspr1 rabbit monoclonal antibody for WB and ELISA, targeting human CNTNAP1. Validated for studying paranodal junctions and myelination in brain tissue. 156 kDa.

[Learn More](#)

Parameter	Value
Product Name	Caspr1 Rabbit mAb
Remarks/Alias	CHN3; P190; CASPR; NRXN4; CNTNAP; Caspr1
Species	Human
GeneID (Human)	8506
GeneID	8506
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 700-800 of human Caspr1 (P78357).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P78357
Protein Weight	156kDa
Shipping	Ice bag