


Anti-Caspr2/Cntnap2 Monoclonal Antibody For Wb, Elisa - Q9Uhc6

Item Number: CM0002500



Introduction

High-quality anti-Caspr2/CNTNAP2 rabbit monoclonal antibody validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying axo-glial junction organization and neurological diseases.

[Learn More](#)

Parameter	Value
Product Name	Caspr2/CNTNAP2 Rabbit mAb
Remarks/Alias	CDFE; NRXN4; AUTS15; CASPR2; PTHSL1; Caspr2/CNTNAP2
Species	Human
GeneID (Human)	26047
GeneID	26047
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1232-1331 of human Caspr2/CNTNAP2 (Q9UHC6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UHC6
Protein Weight	148kDa
Shipping	Ice bag