

Anti-Casq2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O14958

Item Number: CM0020101



Introduction

Rabbit polyclonal antibody against human CASQ2 (Calsequestrin-2), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cardiac excitation-contraction coupling and calcium signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CASQ2 Rabbit pAb
Remarks/Alias	PDIB2; CASQ2
Species	Human
GeneID (Human)	845
GeneID	845
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-399 of human CASQ2 (NP_001223.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14958
Protein Weight	46kDa
Shipping	Ice bag