

Anti-Casr (Extracellular Calcium-Sensing Receptor) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P41180

Item Number: CM0007936



Introduction

Rabbit monoclonal CaSR antibody validated for WB, IF-P, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for calcium signaling and homeostasis research. UniProt P41180.

[Learn More](#)

Parameter	Specification
Product Name	CaSR Rabbit mAb
Remarks/Alias	CAR; FHH; FIH; HHC; EIG8; HHC1; NSHPT; PCAR1; hCasR; GPRC2A; HYPOC1
Species	Human
Gene ID (Human)	846
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 211-310 of human CaSR (NP_000379.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41180
Protein Weight	121 kDa
Shipping	Ice bag