

Anti-Cast Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P20810

Item Number: CM0023328

Introduction

Rabbit polyclonal antibody targeting human CAST (Calpastatin) for Western blot, IHC-P, and ELISA applications. Recombinant immunogen, validated reactivity with human samples. Supports calpain inhibition and muscle protein research.

[Learn More](#)

Parameter	Value
Product Name	CAST Rabbit pAb
Remarks/Alias	BS-17; PLACK; CAST
Species	Human
GeneID (Human)	831
GeneID	831
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human CAST (NP_001177371.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P20810
Protein Weight	77kDa
Shipping	Ice bag