


Anti-Catalase Monoclonal Antibody For Wb, Ihc-P, Elisa - P04040

Item Number: CM0002159



Introduction

Catalase Rabbit Monoclonal Antibody (P04040) validated for WB, IHC-P, ELISA in Human, Mouse, Rat. Targets peroxisomal catalase, an antioxidant enzyme degrading hydrogen peroxide. Ideal for oxidative stress and cell growth research.

[Learn More](#)

Parameter	Value
Product Name	Catalase Rabbit mAb
Remarks/Alias	CAT; catalase; Catalase
Species	Human
Gene ID (Human)	847
Gene ID	847
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Catalase (P04040).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04040
Protein Weight	60kDa
Shipping	Ice bag