


Anti-Catalase Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04040

Item Number: CM0017375



Introduction

High-quality rabbit polyclonal antibody targeting human catalase (CAT), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. The antibody is produced in rabbit and shipped on ice. Suitable for oxidative stress and peroxisomal research.

[Learn More](#)

Parameter	Value
Product Name	Catalase Rabbit pAb
Remarks/Alias	CAT; catalase; MGC138422; MGC138424; Catalase
Species	Human
GenelD (human)	847
GenelD	847
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human Catalase (NP_001743.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04040
Protein Weight	60kDa
Shipping	Ice bag