


Anti-Catalase Rabbit Polyclonal Antibody For Wb And Elisa - P04040

Item Number: CM0013189



Introduction

Rabbit polyclonal antibody to Catalase (CAT), validated for WB and ELISA in human, mouse, and rat samples. Ideal for oxidative stress and cell proliferation research.

[Learn More](#)

Parameter	Specification
Product Name	Catalase Rabbit pAb
Remarks/Alias	CAT; catalase; MGC138422; MGC138424; Catalase
Species	Human
GeneID (human)	847
GeneID	847
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human Catalase (NP_001743.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04040
Protein Weight	60kDa
Shipping	Ice bag