

Anti-Catenin Alpha-2 (Ctnna2) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26232

Item Number: CM0014855



Introduction

Anti-Catenin alpha-2 (CTNNA2) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Recognizes amino acids 167-290. For neural adhesion and development studies.

[Learn More](#)

Parameter	Value
Product Name	CTNNA2 Rabbit pAb
Remarks/Alias	CAPR; CTNR; CAP-R; CT114; CDCBM9; CTNNA2
Species	Human
GeneID (human)	1496
GeneID	1496
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 167-290 of human CTNNA2 (NP_001269526.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26232
Protein Weight	105kDa
Shipping	Ice bag