


Anti-Cathepsin B Monoclonal Antibody For Wb, If/Icc, Elisa - P07858

Item Number: CM0004883



Introduction

Rabbit monoclonal antibody targeting human Cathepsin B (CTSB, APPS), validated for WB, IF/ICC, and ELISA. Detects the ~38 kDa thiol protease involved in protein turnover, tumor invasion, and lysosomal degradation.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin B Rabbit mAb
Remarks/Alias	KWE; APPS; CPSB; RECEUP; Cathepsin B
Species	Human
GeneID (Human)	1508
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 8-339 of human Cathepsin B.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P07858
Protein Weight	38 kDa
Shipping	Ice bag