


Anti-Cathepsin B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07858

Item Number: CM0012877



Introduction

Rabbit polyclonal antibody targeting human Cathepsin B, suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for research in protease biology and cancer metastasis.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin B Rabbit pAb
Remarks/Alias	KWE; APPS; CPSB; RECEUP; Cathepsin B
Species	Human
GeneID (Human)	1508
GeneID	1508
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-339 of human Cathepsin B (NP_001899.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07858
Protein Weight	38kDa
Shipping	Ice bag